

抗ヒト PSA モノクローナル抗体

品 番	K R - 0 3 3
クローン	N o . 6 8
型	I g G 2 a, κ
由 来	B A L B / C マウス腹水
精 製 法	硫酸塩析, r - プロテインA
溶 液	P B S (0.09%アジ化ナトリウム含む)
数 量	1 0 0 μ g
保 存 法	4 ~ 8 °C 凍結融解は出来るだけ避けて下さい。
備 考	抗ヒト PSA 抗体 (No.68) は complex-PSA を測定する際のサンドイッチ ELISA 法の固相に適しています。なお、中和抗体として抗ヒト PSA 抗体 (2E2) および二次抗体としては抗ヒト PSA 抗体 (4D10) が適しています。
参考文献	J Immunol Methods. 2000 Feb 3;234(1-2):99-106.



人と科学のステキな未来へ
コスモ・バイオ株式会社

— 商品に関するお問い合わせ —

TEL: 03-5632-9610 FAX: 03-5632-9619

本社所在地 〒135-0016 東京都江東区東葛 2-2-20 東葛駅前ビル

**MOUSE ANTI PSA
MONOCLONAL ANTIBODY**

Catalog No.	LNМ-KR-033
CLONE:	No.68
QUANTITY:	0.1mg
SPECIFICITY:	Human Prostate specific antigen (PSA) of both free form and complex from.
IMMUNOGEN:	Purified antigen from human seminal fluid
Type/Subtype:	IgG2a κ
APPLICATIONS:	Complex-PSA sandwich ELISA
FORMAT:	Purified immunoglobulin from BALB/C Mouse Ascites (Affinity chromatography)
PRESENTATION:	Liquid in 15mM Phosphate buffer containing 0.15M NaCl, 0.09% sodium azide, pH 7.2
STORAGE/HANDLING:	2 to 8°C, Avoid repeated freeze/thaw cycle.
REFERENCE:	J Immunol Methods. 2000 Feb 3;234(1-2):99-106.

For research use only; not for use as a diagnostic



COSMO BIO CO., LTD.
Inspiration for Life Science

TOYO EKIMAE BLDG. 2-20, TOYO 2-CHOME, KOTO-KU, TOKYO 135-0016 JAPAN
TEL: (81)3-5632-9617 / FAX: (81)3-5632-9618 / e-mail: export@cosmobio.co.jp / URL: www.cosmobio.com