

抗ヒトトランスフェリンレセプターモノクローナル抗体

品 番	K R - 0 4 4
クローン	N o . 5 7
型	I g G 1 , κ
由 来	B A L B / C マウス腹水
精 製 法	硫酸塩析, r - プロテインA
溶 液	P B S (0.09%アジ化ナトリウム含む)
数 量	1 0 0 μ g
保 存 法	4 ~ 8 °C 凍結融解は出来るだけ避けて下さい。
備 考	抗ヒトトランスフェリンレセプター抗体 (No.57) はサンドイッチ ELISA 法の固相に適しています。なお、二次抗体としては抗ヒトトランスフェリンレセプター抗体 (No.143) が適しています。



人と科学のステキな未来へ

コスモ・バイオ株式会社

— 商品に関するお問い合わせ —

TEL: 03-5632-9610 FAX: 03-5632-9619

本社所在地 〒135-0016 東京都江東区東陽 2-2-20 東陽駅前ビル

**MOUSE ANTI TRANSFERRIN RECEPTOR
MONOCLONAL ANTIBODY**

Catalog No.	LNМ-KR-044
CLONE:	No.57
QUANTITY:	0.1mg
SPECIFICITY:	Human Transferrin receptor
IMMUNOGEN:	Human Transferrin receptor
Type/Subtype:	IgG1 κ
APPLICATIONS:	Sandwich ELISA (This antibody can be used as a capture antibody in a human Transferrin receptor ELISA in combination with biotinlated, monoclonal antibody No.143.)
FORMAT:	Purified immunoglobulin from BALB/C Mouse Ascites (Affinity chromatography)
PRESENTATION:	Liquid in 15mM Phosphate buffer containing 0.15M NaCl, 0.09% sodium azide, pH 7.2
STORAGE/HANDLING:	2 to 8°C, Avoid repeated freeze/thaw cycle.

For research use only; not for use as a diagnostic.



COSMO BIO Co., LTD.
Inspiration for Life Science

TOYO EKIMAE BLDG. 2-20, TOYO 2-CHOME, KOTO-KU, TOKYO 135-0016 JAPAN
TEL: (81)3-5632-9617 / FAX: (81)3-5632-9618 / e-mail: export@cosmobio.co.jp / URL: www.cosmobio.com