

抗ヒトTSHモノクローナル抗体

品 番	K R - 0 4 7
クローン	9 G 1 1
型	I g G 1, κ
由 来	B A L B / C マウス腹水
精 製 法	硫酸塩析, r - プロテインA
溶 液	P B S (0.09%アジ化ナトリウム含む)
数 量	1 0 0 μ g
保 存 法	4 ~ 8 °C 凍結融解は出来るだけ避けて下さい。
備 考	抗ヒト TSH 抗体 (9G11) はサンドイッチ ELISA 法の二次抗体に適しています。なお、固相としては抗ヒト TSH 抗体 (3F10) が適しています。



人と科学のステキな未来へ
コスモ・バイオ株式会社

— 商品に関するお問い合わせ —

TEL: 03-5632-9610 FAX: 03-5632-9619

本社所在地 〒135-0016 東京都江東区東葛 2-2-20 東葛駅前ビル

**MOUSE ANTI TSH
MONOCLONAL ANTIBODY**

Catalog No.	LNМ-KR-047
CLONE:	9G11
QUANTITY:	0.1mg
SPECIFICITY:	Human Thyroid stimulating hormone (β -chain specific) No cross reactivity to HCG
IMMUNOGEN:	Isolated from human pituitary glands
Type/Subtype:	IgG1 κ
APPLICATIONS:	Sandwich ELISA (This antibody can be used as a detection antibody in a human TSH ELISA combination with the monoclonal capture antibody 3F10.)
FORMAT:	Purified immunoglobulin from BALB/C Mouse Ascites (Affinity chromatography)
PRESENTATION:	Liquid in 15mM Phosphate buffer containing 0.15M NaCl, 0.09% sodium azide, pH 7.2
STORAGE/HANDLING:	2 to 8°C, Avoid repeated freeze/thaw cycle.

For research use only; not for use as a diagnostic.



COSMO BIO CO., LTD.
Inspiration for Life Science

TOYO EKIMAE BLDG. 2-20, TOYO 2-CHOME, KOTO-KU, TOKYO 135-0016 JAPAN
TEL: (81)3-5632-9617 / FAX: (81)3-5632-9618 / e-mail: export@cosmobio.co.jp / URL: www.cosmobio.com