



COSMO BIO CO., LTD.
Inspiration for Life Science

MOUSE ANTI ADAMTS13 MONOCLONAL ANTIBODY

Catalog No.	LNM-KR-052
Clone	C7
Concentration	1 mg/ml
Volume	100 UL
Subclass	IgG1k
Immunogen	recombinant ADAMTS13
Specificity	ADAMTS13 [Tsp7-8 domain]
Applications	Suitable for immunoblotting (non-reducing condition), Immunohistochemistry, ELISA (detection antibody)
Format	Purified immunoglobulin from BALB/C Mouse Ascites Ammonium sulfate fractionation Affinity chromatography
Formulation	Liquid in 15mM Phosphate buffer containing 0.15M NaCl, 0.09% sodium azide, pH 7.2
Storage/Handling	2 to 8°C, Maintain at -70°C for up to 12 months. Avoid repeated freeze/thaw cycle.
Reference	Uemura M et al., Blood., 106, 922-924 (2005)

For research use only; not for use as a diagnostic



COSMO BIO CO., LTD.
Inspiration for Life Science

TOYO EKIMAE BLDG, 2-20, TOYO 2CHOME, KOTO-KU, TOKYO 135-0016, JAPAN
TEL: +81-3-5632-9617 FAX: +81-3-5632-9618 e-mail: export@cosmobio.co.jp

www.cosmobio.co.jp

抗 ADAMTS13 モノクローナル抗体

品 番	K R - 0 5 2
クローン	C 7
型	I g G 1, κ
免疫抗原	recombinant ADAMTSA13
特異性	ADAMTS13 (Tsp7-8 domain)
由 来	B A L B / C マウス腹水
精 製 法	硫酸塩析, r - プロテイン A
溶 液	P B S (0.09% アジ化ナトリウム含む)
数 量	1 0 0 μ g
保 存 法	4 ~ 8 °C 長期保存の場合は、- 7 0 °C 以下にて凍結保存して下さい。 凍結融解は出来るだけ避けて下さい。
備 考	WB (非還元条件)、IH、ELISA (二次抗体)
Reference	Uemura M et al., Blood., 106., 922-924.(2005)



コスモ・バイオ株式会社