

## 抗ヒト Ca125 モノクローナル抗体

|       |                                                                                         |
|-------|-----------------------------------------------------------------------------------------|
| 品 番   | K R - 0 6 0                                                                             |
| クローン  | N o . 5 2                                                                               |
| 型     | I g G 1 , $\kappa$                                                                      |
| 由 来   | B A L B / C マウス腹水                                                                       |
| 精 製 法 | 硫酸塩析, r - プロテインA                                                                        |
| 溶 液   | P B S (0.09%アジ化ナトリウム含む)                                                                 |
| 数 量   | 1 0 0 $\mu$ g                                                                           |
| 保 存 法 | 4 ~ 8 °C<br>凍結融解は出来るだけ避けて下さい。                                                           |
| 備 考   | 抗ヒト Ca125 抗体 (No.52) はサンドイッチ ELISA 法の固相に適しています。なお、二次抗体としては抗ヒト Ca125 抗体 (No.35) が適しています。 |



人と科学のステキな未来へ  
**コスモ・バイオ株式会社**

— 商品に関するお問い合わせ —

TEL: 03-5632-9610 FAX: 03-5632-9619

本社所在地 〒135-0016 東京都江東区東葛 2-2-20 東葛駅前ビル

**MOUSE ANTI Ca125  
MONOCLONAL ANTIBODY**

|                          |                                                                                                                                                      |
|--------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Catalog No.</b>       | LNМ-KR-060                                                                                                                                           |
| <b>CLONE:</b>            | No.52                                                                                                                                                |
| <b>QUANTITY:</b>         | 0.1mg                                                                                                                                                |
| <b>SPECIFICITY:</b>      | Human Ca125                                                                                                                                          |
| <b>IMMUNOGEN:</b>        | Human Ca125                                                                                                                                          |
| <b>Type/Subtype:</b>     | IgG1 $\kappa$                                                                                                                                        |
| <b>APPLICATIONS:</b>     | Sandwich ELISA ( This antibody can be used as a capture antibody in a human Ca125 ELISA in combination with biotinlated, monoclonal antibody No.35.) |
| <b>FORMAT:</b>           | Purified immunoglobulin from BALB/C Mouse Ascites<br>(Affinity chromatography)                                                                       |
| <b>PRESENTATION:</b>     | Liquid in 15mM Phosphate buffer containing 0.15M NaCl,<br>0.09% sodium azide, pH 7.2                                                                 |
| <b>STORAGE/HANDLING:</b> | 2 to 8°C, Avoid repeated freeze/thaw cycle.                                                                                                          |

**For research use only; not for use as a diagnostic.**



**COSMO BIO CO., LTD.**  
Inspiration for Life Science

TOYO EKIMAE BLDG. 2-20, TOYO 2-CHOME, KOTO-KU, TOKYO 135-0016 JAPAN  
TEL: (81)3-5632-9617 / FAX: (81)3-5632-9618 / e-mail: export@cosmobio.co.jp / URL: www.cosmobio.com