

## Anti (human) Fibrinectin

### 1. Description

<b>Host animal</b>	Rabbit
<b>Source (Volume)</b>	Whole plasma(100ul)
<b>Titer</b>	According to the ELISA assay, results are positive for dilutions up to 100,000 fold.
<b>Source of antigen</b>	human plasma
<b>Cross reactivity</b>	Cross react with guinea pig, bovine and pig.*
<b>Characteristic</b>	Does not cross react with type I, II, III, IV, V collagen and laminin.
<b>Standard dilution</b>	More than 1:1,000 dilution using a chemiluminescent method .

\* Cross reaction will differ between species. The above dilution is only a recommendation and the optimum concentration may differ for each case.

### 2. Storage

Store below -20°C ( below -70°C for prolonged storage).  
After thawing, store in small aliquots in sealable vials and store below -70°C. To prevent degradation from repeated thawing, store the antiserum between 0 to 4°C after second thawing.

### 3. Stability

Stable for three years at -70°C.  
This product does not contain preservatives such as NaN3.

*For research use only; not for use as a diagnostic.*

**COSMO BIO CO.,LTD**

TOYO EKIMAE BLDG, 2-20,TOYO 2CHOME,KOTO-KU,TOKYO 135-0016,JAPAN  
TEL: +81-3-5632-9605 FAX: +81-3-5632-9614 e-mail: export@cosmobio.co.jp

[www.cosmobio.co.jp](http://www.cosmobio.co.jp)