

Anti (human) Fibrinectin

1. Description

Host animal	Rabbit
Source (Volume)	Whole plasma(100ul)
Titer	According to the ELISA assay, results are positive for dilutions up to 100,000 fold.
Source of antigen	human plasma
Cross reactivity	Cross react with guinea pig, bovine and pig.*
Characteristic	Does not cross react with type I, II, III, IV, V collagen and laminin.
Standard dilution	More than 1:1,000 dilution using a chemiluminescent method .

* Cross reaction will differ between species. The above dilution is only a recommendation and the optimum concentration may differ for each case.

2. Storage

Store below -20°C (below -70°C for prolonged storage).
After thawing, store in small aliquots in sealable vials and store below -70°C. To prevent degradation from repeated thawing, store the antiserum between 0 to 4°C after second thawing.

3. Stability

Stable for three years at -70°C.
This product does not contain preservatives such as NaN3.

For research use only; not for use as a diagnostic.

COSMO BIO CO.,LTD

TOYO EKIMAE BLDG, 2-20,TOYO 2CHOME,KOTO-KU,TOKYO 135-0016,JAPAN
TEL: +81-3-5632-9605 FAX: +81-3-5632-9614 e-mail: export@cosmobio.co.jp

www.cosmobio.co.jp