



COSMO BIO Co., LTD.

Inspiration for Life Science

Catalog No. LSL-LB-1027-EX

Anti (bovine) Fibronectin

1. Description

| | |
|--------------------------|--|
| Host animal | Rabbit |
| Source (Volume) | Whole serum (100ul) |
| Titer | According to the ELISA assay, results are positive for dilutions up to 300,000 fold. |
| Source of antigen | bovine plasma |
| Cross reactivity | Cross react with bovine, goat, guinea pig, rat, mouse and human.* |
| Characteristic | Does not cross react with bovine type I, II, III, IV, V, VI collagen and mouse laminin. |
| Standard dilution | More than 1:500 dilution using a chemiluminescent method and more than 1:1,000 dilution using EIA. |

* Cross reaction will differ between species. The above dilution is only a recommendation and the optimum concentration may differ for each case.

2. Storage

Store below -20°C (below -70°C for prolonged storage).
After thawing, store in small aliquots in sealable vials and store below -70°C. To prevent degradation from repeated thawing, store the antiserum between 0 to 4°C after second thawing.

3. Stability

Stable for three years at -70°C.
This product does not contain preservatives such as NaN3.

For research use only; not for use as a diagnostic.



COSMO BIO Co., LTD.
Inspiration for Life Science

TOYO EKIMAE BLDG, 2-20, TOYO 2CHOME, KOTO-KU, TOKYO 135-0016, JAPAN
TEL: +81-3-5632-9617 FAX: +81-3-5632-9618 e-mail: export@cosmobio.co.jp

www.cosmobio.co.jp