

Anti (human) Type V Collagen

1. Description

Host animal	Rabbit
Source (Volume)	Whole serum (100ul)
Titer	According to the ELISA assay, results are positive for dilutions up to 50,000 fold.
Source of antigen	human placenta
Cross reactivity	Cross react with human.
Characteristic	Does not cross react with type I, II, III, IV, VI collagen.
Standard dilution	More than 1:300 dilution (recommended : 1/200 ~ 1/1,000)

using
a chemiluminescent method.

2. Storage

Store below -20°C (below -70°C for prolonged storage).
After thawing, store in small aliquots in sealable vials and store below -70°C. To prevent degradation from repeated thawing, store the antiserum between 0 to 4°C after second thawing.

3. Stability

Stable for three years at -70°C.
This product does not contain preservatives such as NaN₃.

For research use only; not for use as a diagnostic.

COSMO BIO CO.,LTD

TOYO EKIMAE BLDG, 2-20, TOYO 2CHOME, KOTO-KU, TOKYO 135-0016, JAPAN
TEL: +81-3-5632-9605 FAX: +81-3-5632-9614 e-mail: export@cosmobio.co.jp

www.cosmobio.co.jp