



COSMO BIO Co., LTD.

Inspiration for Life Science

Catalog No. LSL-LB-2007-EX

## Anti (bovine) Vitronectin

### 1. Description

<b>Host animal</b>	Rabbit
<b>Source (Volume)</b>	Whole serum (100ul)
<b>Titer</b>	According to the ELISA assay, results are positive for dilutions up to 500,000 fold.
<b>Source of antigen</b>	Bovine plasma
<b>Cross reactivity</b>	Cross react with bovine.
<b>Characteristic</b>	Does not cross react with Type I, II, III, IV, V, VI collagen, fibronectin, laminin, elastin and serum albumin. Appears as positive band by Western blotting.
<b>Standard dilution</b>	More than 1:500 dilution (recommended: 1/500~1/2,000)* using a chemiluminescent method and more than 1:2,000 dilution using Western Blotting.

\* Cross reaction will differ between species. The above dilution is only a recommendation and the optimum concentration may differ for each case.

### 2. Storage

Store below -20°C ( below -70°C for prolonged storage).  
After thawing, store in small aliquots in sealable vials and store below -70°C. To prevent degradation from repeated thawing, store the antiserum between 0 to 4°C after second thawing.

### 3. Stability

Stable for three years at -70°C.  
This product does not contain preservatives such as NaN<sub>3</sub>.

*For research use only; not for use as a diagnostic.*



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO EKIMAE BLDG, 2-20, TOYO 2CHOME, KOTO-KU, TOKYO 135-0016, JAPAN  
TEL: +81-3-5632-9617 FAX: +81-3-5632-9618 e-mail: export@cosmobio.co.jp

[www.cosmobio.co.jp](http://www.cosmobio.co.jp)