



Anti (chMAP) Myelin protein zero (P0)

Form	Whole serum
Host	Rabbit
Antigen	Multiple Antigen Peptide*
Application notes	ELISA/ IHC/ WB Recommended use and dilutions According to the ELISA assay, results are positive for dilutions up to 200,000 fold against linear synthetic peptides For Western Blotting, dilute > 1 :1,000 For Immunofluorescence, dilute> 1 :500 (1/200 – 1/2000) Optimal diluting should be determined by the user.
Cross Reactivity	HU/ MS/ RAT/ CHK
Volume	100 ul
Storage	Store below -40°C (Store below -80°C for long term storage) After thawing, Avoid freeze-thaw.
Stability	For 3 years at -80C Preservatives such as NaN ₃ are not included

*The cross reactivity is different depending on the species. Please confirm the dilution ratio.

