



Anti Mite Extract

1. Description

Host animal	Rabbit
Source (Volume)	Whole serum (100ul)
Source of antigen	Mite Crude Extract (#LG-5119)
Characteristic	Positive to 25,000 times dilution Mite Extract-Df (#LG-5339) and to 50,000 times dilution Mite Extract-Dp (#LG-5449) and to 3,000 times dilution Mite Extract-Tp (#LG-5559) .

2. Storage

Store below -20°C (below -70°C for prolonged storage).
After thawing, store in small aliquots in sealable vials and store below -70°C. To prevent degradation from repeated thawing, store the antiserum between 0°C to 4°C after second thawing.

3. Stability

Stable for three years at -70°C.
This product does not contain preservatives such as NaN₃.

For research use only; not for use as a diagnostic.

