



Anti TNP-Ascaris

1. Description

Host animal	Rat
Source (Volume)	Whole serum (100ul)
Source of antigen	2,4,6-trinitrophenyl Ascaris Extract (#LG-1109)
Standard dilution	Positive more than 1:500 dilution using 48h rat homologous passive cutaneous anaphylaxis (PCA). 1999, 06, 26 Official approval

2. Storage

Store below -20°C (below -70°C for prolonged storage).
After thawing, store in small aliquots in sealable vials and store below -70°C. To prevent degradation from repeated thawing, store the antiserum between 0 to 4°C after second thawing.

3. Stability

Stable for three years at -70°C.
This product does not contain preservatives such as NaN₃.

For research use only; not for use as a diagnostic.

