

POLYCLONAL ANTIBODY

For research use only. Not for clinical diagnosis.

Catalog No. LSL-LB-9009

Anti DNP-Ascaris

Product type	Primary antibodies
Host	Rat
Source	Serum
Form	Liquid
Volume	This product does not contain preservatives such as NaN ₃ .
Concentration	100μl
Specificity	DNP-Ascaris
Antigen	2,4,-dinitrophenyl Ascaris Extract (#LG-0009)
Isotype	

Application notes ELISA, WB, IF, IHC
Recommended use
Using 48h rat homologous passive cutaneous anaphylaxis (PCA).

Immunohistochemistry on Frozen or Paraffin section

Recommended dilutions

Optimal dilutions/concentrations should be determined by the end user.

Staining Pattern

Cross reactivity Cross reaction will differ between species.

Storage Store below -20°C (below -70°C for prolonged storage). After thawing, store in small aliquots in sealable vials and store below -70°C. To prevent degradation from repeated thawing, store the antiserum between 0 to 4°C after second thawing.

Stable for three years at -70°C

References

For research use only. Not for clinical diagnosis.



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619

Anti DNP-Ascaris

I. 内容

免疫動物	ラット
性状・包装サイズ	全血清・100 μ l
抗原由来	2,4-dinitrophenylated Ascaris Extract (LG-0009)
標準希釈率	48 hr ラット homologous passive cutaneous anaphylaxis (PCA) で1:512以上

II. 保存上の注意

−20℃以下（長期間の場合は−70℃以下）で凍結して下さい。
解凍後は密栓のできる小型容器に研究の規模に応じて少量ずつ分注し、−70℃以下で保存して下さい。
凍結融解の繰り返しによる力価の低下を避ける為、再解凍後の抗血清は0〜4℃に保ち操作・保存して下さい。

III. 安定性

−70℃で3年間安定。
但し、NaN₃等の防腐剤は入っていません。

製造元

総発売元

LSL 株式会社 エル・エス・エル

コスモ・バイオ株式会社