



Anti DNP-Ascaris

1. Description

| | |
|--------------------------|---|
| Host animal | Rat |
| Source (Volume) | Whole serum (100ul) |
| Source of antigen | 2,4,-dinitrophenyl Ascaris Extract (#LG-0009) |
| Standard dilution | using 48h rat homologous passive cutaneous anaphylaxis (PCA). |

2. Storage

Store below -20°C (below -70°C for prolonged storage).
After thawing, store in small aliquots in sealable vials and store below -70°C. To prevent degradation from repeated thawing, store the antiserum between 0 to 4°C after second thawing.

3. Stability

Stable for three years at -70°C.
This product does not contain preservatives such as NaN₃.

For research use only; not for use as a diagnostic.

