



DNP-OA (Ovalbumin)

| | |
|----------------------|---|
| Form | Lyophilized protein powder 2,4-dinitrophenylated Chicken Egg Albumin |
| Volume | Protein: 30mg Conjugation ratio: 3.6 DNP groups/molecule (Lot #802051) (As average of molecular weight: 45,000) |
| Concentration | Add 2.0mL of distilled water. It will become 5mM Borate buffered saline [BBS] (pH 8.0), 0.9% NaCl solution. |
| Storage | Store below 4°C (below -20°C for prolonged storage or after dissolved) |

Allergic antigen

| | |
|---------|--------------------------------------|
| LG-0009 | DNP-Ascaris |
| LG-0533 | DNP-Mite-Df |
| LG-0544 | DNP-Mite-Dp |
| LG-0528 | DNP-Cedar Pollen-Cj |
| LG-0017 | DNP-BSA |
| LG-3017 | DNP-BSA (30) |
| LG-0024 | DNP-OA (Ovalbumin) |
| LG-0047 | DNP-Casein |
| LG-0067 | DNP-LG (BETA-Lactoglobulin) |
| LG-0089 | DNP-KLH |
| LG-0317 | BPO-BSA |
| LG-0577 | BPO-BGG |
| LG-1109 | TNP-Ascaris |
| LG-1117 | TNP-BSA |
| LG-5009 | Ascaris Extract |
| LG-5339 | Mite Extract-Df |
| LG-5449 | Mite Extract-Dp |
| LG-5559 | Mite Extract-Tp |
| LG-5229 | Cedar Pollen Extract |
| LG-5280 | Cedar Pollen Extract-Cj |
| LG-5779 | Hinoki Cypress, Pollen Crude Extract |
| LG-5889 | Ragweed Pollen Extract |

Adjuvant

| | |
|---------|------|
| LG-6000 | ALUM |
|---------|------|

For research use only. Not for clinical diagnosis.



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619

DNP-0A (Ovalbumin)

I. 内容

lot No. 802051

| | |
|-----|---|
| 性 状 | 2,4-dinitrophenylated Chicken Egg Albumin 凍結乾燥粉末 |
| 規 格 | 蛋白質量：30mg DNP基数：1分子あたり3.6個（平均分子量45,000 としての計算値） 2.0mlの精製水を加えることで5mM硼酸緩衝溶液 （pH8.0）0.9%NaCl溶液になります。 |

II. 保存上の注意

4℃以下（長期間又は溶液の場合は-20℃以下）で保存して下さい。

製造元

総発売元

 株式会社 エル・エス・エル

コスモ・バイオ株式会社

新製品案内

| | |
|--------------------------|-----------------------|
| ☆Mite Extract | カタログNo. LG-5119 (LSL) |
| ☆Mite Extract-Df | カタログNo. LG-5339 (LSL) |
| ☆Mite Extract-Dp | カタログNo. LG-5449 (LSL) |
| ☆Mite Extract-Tp | カタログNo. LG-5559 (LSL) |
| ☆Cedar Pollen Extract | カタログNo. LG-5229 (LSL) |
| ★Cedar Pollen Extract-Cj | カタログNo. LG-5280 (LSL) |
| ★FTPP (フォルマリン処理百日咳死菌) | カタログNo. LG-7000 (LSL) |