



Mite Extract-Dp

1. Description

Freezed Dry powder (Mite-Dp Crude Extract)

This product extracted from House Dust Mites, Dermatophagoides pteronyssinus (Trouessari)*.

Protein : 10mg

Add 2.0ml of distilled water. It will become 5mM Borate buffered saline (pH8.0), 0.9% NaCl [BBS] solution.

* Mites were homogenized with 10-fold 50mM Phosphate buffer (pH7.2) in an ice bath, and the homogenate were stirred at 4°C for 24 hours, and centrifuged at dialyzed 4 times by BBS.

2. Storage

Store below 4°C (below -20°C for prolonged storage).

For research use only; not for use as a diagnostic.

