

MAB Institute, Inc.
82 Inaokacho, Yokosuka-shi,
Kanagawa 238-8580, Japan
TEL, FAX: +81 46-822-9328

Anti-acetyl Histone H3 (Lys27), mouse monoclonal antibody

Catalog No. MABI0009-20, MABI0009-100

Lot No. 09001

Product Description

Clone No.	MABI0309 (CMA309)
Host	Mouse
Source	Culture supernatant of serum free medium
Isotype	IgG1
Size	100 μ l
Concentration	1 mg/ml
Antigen	19 amino acid residues around the Lys27 of human Histone H3.1
Buffer	PBS+0.05% sodium azide
Storage	Store below -20°C
Expiration Date	Mar 2011

Application notes

Immunoprecipitation	1-5 μ g/5 μ l Sepharose
Immunostaining	0.5-1 μ g/ml
Immunoblot	0.5-1 μ g/ml

Attention

Mouse IgG has weak combination with protein A and protein G. Please use the beads which connect anti-mouse IgG antibodies such as anti-mouse IgG sepharose for immunoprecipitation.

Reference

Kimura H, Hayashi-Takanaka Y, Goto Y, Takizawa N, Nozaki N. The organization of histone H3 modifications as revealed by a panel of specific monoclonal antibodies. Cell Struct Funct. 2008;33(1):61-73.

MABI 0309

