

## **Monoclonal Antibody**

### **to Human Chorionic Gonadotropin(MCG-040)**

<b>CATALOG NUMBER:</b>	NBT-MCG-040
<b>QUANTITY:</b>	1mg
<b>CLONE:</b>	CG-040
<b>SOURCE OF ANTIGEN:</b>	Human Chorionic Gonadotropin
<b>SOURCE:</b>	BALB/C Mouse implantation ascites
<b>FORMAT:</b>	Purified Protein
<b>CONCENTRATION:</b>	10mg/ml
<b>ISOTYPE:</b>	IgG1k
<b>PURIFY:</b>	Protein A
<b>PURITY:</b>	1 band (IgG1k) Cellulose acetate Electrophoresis
<b>PRESENTATION:</b>	Liquid in 10mM phosphate buffer with 0.15M NaCl , pH7.4, 0.1% Sodium Azide
<b>STORAGE/HANDLING:</b>	Maintain at 4 °C-8 °C (-20 °C for long term storage) in undiluted aliquots. Avoid repeated freeze/thaw cycles

*For research use only; not for use as a diagnostic.*