

Monoclonal Antibody
to Human CRP[C-Reactive Protein](MCR-019)

CATALOG NUMBER:	NBT-MCR-019
QUANTITY:	1mg
CLONE:	CRB-019
SOURCE OF ANTIGEN:	Human CRP[C-Reactive Protein]
SOURCE:	BALB/C Mouse implantation ascites
FORMAT:	Purified Protein
CONCENTRATION:	10mg/ml
ISOTYPE:	IgG2b
PURIFY:	Protein A
PURITY:	1 band (IgG1k) Cellulose acetate Electrophoresis
PRESENTATION:	Liquid in 10mM phosphate buffer with 0.15M NaCl , pH7.4, 0.1% Sodium Azide
STORAGE/HANDLING:	Maintain at 4 °C-8 °C (-20 °C for long term storage) in undiluted aliquots. Avoid repeated freeze/thaw cycles.

For research use only; not for use as a diagnostic.