



COSMO BIO CO., LTD.
Inspiration for Life Science

Monoclonal Antibody

to Human CRP[C-Reactive Protein](MCR-020)

CATALOG NUMBER:	NBT-MCR-020
QUANTITY:	1mg
CLONE:	CRB-020
SOURCE OF ANTIGEN:	Human CRP[C-Reactive Protein]
SOURCE:	BALB/C Mouse implantation ascites
FORMAT:	Purified Protein
CONCENTRATION:	10mg/ml
ISOTYPE:	IgG2b
PURIFY:	Protein A
PURITY:	1 band (IgG1k) Cellulose acetate Electrophoresis
PRESENTATION:	Liquid in 10mM PBS, pH7.2 0.09% Sodium Azide
STORAGE/HANDLING:	Maintain at 4 °C-8 °C (-20 °C for long term storage) in undiluted aliquots. Avoid repeated freeze/thaw cycles.

For research use only; not for use as a diagnostic.



COSMO BIO CO., LTD.
Inspiration for Life Science

TOYO EKIMAE BLDG, 2-20, TOYO 2CHOME, KOTO-KU, TOKYO 135-0016, JAPAN
TEL: +81-3-5632-9617 FAX: +81-3-5632-9618 e-mail: export@cosmobio.co.jp

www.cosmobio.co.jp

株式会社日本バイオテスト研究所
〒185-0002
東京都国分寺市東戸倉 1-15-3
Tel 042-321-0921
Fax 042-323-7025

略号 NBT
MCR シリーズ
カタログ番号 MCR-017~023
品名 抗ヒト C-反応性蛋白質 モノクローナル抗体

本製品はヒト C-反応性蛋白質を BALB/c マウスに免疫し、その脾細胞とマウス骨髄腫細胞 (P3U1) とをポリエチレングリコール法により融合し得たハイブリドーマを BALB/c マウス腹腔中で増殖させ、その腹水液を精製したマウスモノクローナル抗体です。

抗原 : 精製ヒト C-反応性蛋白質

特性 : 表に示す

カタログ No.	クローン	サブクラス	性 質
MCR-017	CRB-017	IgG1	CRP と凝集反応を示します
MCR-018	CRB-018	IgG1	CRP と凝集反応を示します 反応は EDTA によって阻害されます
MCR-019	CRB-019	IgG1	CRP と凝集反応を示します 反応は EDTA によって阻害されます
MCR-020	CRB-020	IgG2b	CRP と凝集反応を示します
MCR-023	CRB-023	IgG1	CRP と凝集反応を示します 反応は EDTA によって阻害されます

内容 : 10mM PBS pH7.2 (0.09NaN₃) にて 10mg/ml に調整してあります。
保存法 : 4~8℃

承認	作成
年 月 日	2022年6月27日
	