

Monoclonal Antibody to D Dimer (DF-042)

CATALOG NUMBER:	NBT-MDF-042
LOT NUMBER:	010418-1
QUANTITY:	1mg
CLONE:	DF-042
SUBCLASS:	IgG1k
SOURCE:	BALB/C Mouse implantation ascites
FORMAT:	Purified Protein
CONCENTRATION:	10mg/ml
PURIFY:	Protein A
PURITY:	1 band (IgG1k) Cellulose acetate Electrophoresis
PRESENTATION:	Liquid in 0.15M NaCl 10mM PB, pH7.2, 0.1% Sodium Azide
STORAGE/HANDLING:	Maintain at 4 °C-8 °C (-20 °C for long term storage) in undiluted aliquots. Avoid repeated freeze/thaw cycles.

For research use only; not for use as a diagnostic.

COSMO BIO CO.,LTD

TOYO EKIMAE BLDG, 2-20, TOYO 2CHOME, KOTO-KU, TOKYO 135-0016, JAPAN
TEL: +81-3-5632-9605 FAX: +81-3-5632-9614 e-mail: export@cosmobio.co.jp

www.cosmobio.co.jp