

## **Monoclonal Antibody to Human Haptoglobin (MFG-101)**

|                           |   |
|---------------------------|---|
| <b>CATALOG NUMBER:</b>    | NBT-MFG-101   |
| <b>QUANTITY:</b>          | 1mg   |
| <b>SOURCE OF ANTIGEN:</b> | Human haptoglobin   |
| <b>CLONE:</b>             | FG-101  |
| <b>SOURCE OF ANTIGEN:</b> | Human haptoglobin   |
| <b>SOURCE:</b>            | BALB/C Mouse implantation ascites   |
| <b>FORMAT:</b>            | Purified Protein  |
| <b>CONCENTRATION:</b>     | 10mg/ml   |
| <b>ISOTYPE:</b>           | IgG2ak  |
| <b>PURIFY:</b>            | Protein A   |
| <b>PURITY:</b>            | 1 band (IgG1k) Cellulose acetate Electrophoresis  |
| <b>PRESENTATION:</b>      | Liquid in 10mM phosphate buffer with 0.15M NaCl , pH7.4,<br>0.1% Sodium Azide                                     |
| <b>STORAGE/HANDLING:</b>  | maintain at 4 °C-8 °C (-20 °C for long term storage) in undiluted<br>aliquots. Avoid repeated freeze/thaw cycles. |

*For research use only; not for use as a diagnostic.*