



COSMO BIO CO., LTD.
Inspiration for Life Science

Monoclonal Antibody

to Multiubiquitin Chains (FK1)

CATALOG NUMBER:	NBT-MFK001-EX
QUANTITY:	0.5mg
CLONE:	FK1
SOURCE OF ANTIGEN:	crude preparation of polyubiquitinated-lysozyme
SOURCE:	BALB/C Mouse implantation ascites
SPECIFICITY:	Does not react with methyl ubiquitinated protein and free-ubiquitin
FORMAT:	Purified Protein
CONCENTRATION:	10mg/ml
ISOTYPE:	IgM
PURIFY:	DEAE ion exchange column chromatography
PURITY:	1 band (IgM) Cellulose acetate Electrophoresis
PRESENTATION:	Liquid in 10mM phosphate buffer with 0.15M NaCl , pH7.4, 0.1% Sodium Azide
STORAGE/HANDLING:	Maintain at 4 °C-8 °C (-20 °C for long term storage) in undiluted aliquots. Avoid repeated freeze/thaw cycles.
REFERENCE:	(1)M.Fujimoto, H.Sawada, and H. Yokosawa, FEBS Lett.349,173-180(1994) (2)K. Takada, H.Nasu, H,Hibi, Y. Tsukada, K. Ohkawa, M.Fujimuro, H. Sawada, and H. Yokosawa, Eur. J. Biochem. 223, 42-47(1995)

For research use only; not for use as a diagnostic.



COSMO BIO CO., LTD.
Inspiration for Life Science

TOYO EKIMAE BLDG, 2-20, TOYO 2CHOME, KOTO-KU, TOKYO 135-0016, JAPAN
TEL: +81-3-5632-9617 FAX: +81-3-5632-9618 e-mail: export@cosmobio.co.jp

www.cosmobio.co.jp