

Monoclonal Antibody to Human Albumin (MFU-302)

CATALOG NUMBER:	NBT-MFU-302
QUANTITY:	1mg
SOURCE OF ANTIGEN:	Human Albumin
CLONE:	FU-302
SOURCE OF ANTIGEN:	Human Albumin
SOURCE:	BALB/C Mouse implantation ascites
FORMAT:	Purified Protein
CONCENTRATION:	10mg/ml
ISOTYPE:	IgG1k
PURIFY:	DEAE ion exchange column chromatography
PURITY:	1 band (IgG1k) Cellulose acetate Electrophoresis
PRESENTATION:	Liquid in 10mM phosphate buffer with 0.15M NaCl , pH7.4, 0.1% Sodium Azide
STORAGE/HANDLING:	maintain at 4 °C-8 °C (-20 °C for long term storage) in undiluted aliquots. Avoid repeated freeze/thaw cycles.

For research use only; not for use as a diagnostic.