

## **Monoclonal Antibody to Human Albumin (MFU-303)**

<b>CATALOG NUMBER:</b>	NBT-MFU-303
<b>QUANTITY:</b>	1mg
<b>SOURCE OF ANTIGEN:</b>	Human Albumin
<b>CLONE:</b>	FU-303
<b>SOURCE OF ANTIGEN:</b>	Human Albumin
<b>SOURCE:</b>	BALB/C Mouse implantation ascites
<b>FORMAT:</b>	Purified Protein
<b>CONCENTRATION:</b>	10mg/ml
<b>ISOTYPE:</b>	IgG1k
<b>PURIFY:</b>	DEAE ion exchange column chromatography
<b>PURITY:</b>	1 band (IgG1k) Cellulose acetate Electrophoresis
<b>PRESENTATION:</b>	Liquid in 10mM phosphate buffer with 0.15M NaCl , pH7.4, 0.1% Sodium Azide
<b>STORAGE/HANDLING:</b>	maintain at 4 °C-8 °C (-20 °C for long term storage) in undiluted aliquots. Avoid repeated freeze/thaw cycles.

*For research use only; not for use as a diagnostic.*