

Monoclonal Antibody to Human IgG (IG-002)

CATALOG NUMBER:	NBT-MIG-002
QUANTITY:	1mg
CLONE:	IG-002
SOURCE OF ANTIGEN:	Purified Human IgG
SOURCE:	BALB/C Mouse implantation ascites
FORMAT:	Purified Protein
CONCENTRATION:	10mg/ml
ISOTYPE:	IgG1k
PURIFY:	DEAE ion exchange chromatography
PURITY:	1 band (IgG1k) Cellulose acetate Electrophoresis
TITER:	According to the ELISA assay, results are positive for dilutions up to 12,500 fold .
PRESENTATION:	Liquid in 10mM phosphate buffer with 0.15M NaCl , pH7.4, 0.1% Sodium Azide
STORAGE/HANDLING:	Maintain at 4 °C-8 °C (-20 °C for long term storage) in undiluted aliquots. Avoid repeated freeze/thaw cycles.

For research use only; not for use as a diagnostic.

COSMO BIO CO.,LTD

TOYO EKIMAE BLDG, 2-20, TOYO 2CHOME, KOTO-KU, TOKYO 135-0016, JAPAN
TEL: +81-3-5632-9605 FAX: +81-3-5632-9614 e-mail: export@cosmobio.co.jp

www.cosmobio.co.jp