

Monoclonal Antibody to Human Transferrin (TH-003)

| | |
|---------------------------|---|
| CATALOG NUMBER: | NBT-MTH-003 |
| QUANTITY: | 1mg |
| CLONE: | TH-003 |
| SOURCE OF ANTIGEN: | purified holo-transferrin |
| SOURCE: | BALB/C Mouse implantation ascites |
| FORMAT: | purified Protein |
| CONCENTRATION: | 10mg/ml |
| ISOTYPE: | IgG1k |
| PURIFY: | proteinA column chromatography |
| PURITY: | 1 band (IgG1k) Cellulose acetate Electrophoresis |
| PRESENTATION: | Liquid in 10mM phosphate buffer with 0.15M NaCl , pH7.4, 0.1% Sodium Azide |
| STORAGE/HANDLING: | Maintain at 4 °C-8 °C (-20 °C for long term storage) in undiluted aliquots. Avoid repeated freeze/thaw cycles. |

For research use only; not for use as a diagnostic.

COSMO BIO CO.,LTD

TOYO EKIMAE BLDG, 2-20, TOYO 2CHOME, KOTO-KU, TOKYO 135-0016, JAPAN
TEL: +81-3-5632-9605 FAX: +81-3-5632-9614 e-mail: export@cosmobio.co.jp

www.cosmobio.co.jp