

## Monoclonal Antibody to Human Insulin (YO-201)

<b>CATALOG NUMBER:</b>	NBT-MYO-201
<b>QUANTITY:</b>	1mg
<b>CLONE:</b>	YO-201
<b>SOURCE OF ANTIGEN:</b>	Human Insulin
<b>SOURCE:</b>	BALB/C Mouse implantation ascites
<b>FORMAT:</b>	Purified Protein
<b>CONCENTRATION:</b>	10mg/ml
<b>ISOTYPE:</b>	IgG1
<b>PURIFY:</b>	proteinA column chromatography
<b>PURITY:</b>	1 band (IgG1) Cellulose acetate Electrophoresis
<b>TITER:</b>	According to the ELISA assay, results are positive (OD 0.1) for dilutions up to 45,000 fold .
<b>PRESENTATION:</b>	Liquid in 10mM phosphate buffer with 0.15M NaCl , pH7.4, 0.1% Sodium Azide
<b>STORAGE/HANDLING:</b>	Maintain at 4 °C-8 °C (-20 °C for long term storage) in undiluted aliquots. Avoid repeated freeze/thaw cycles.

*For research use only; not for use as a diagnostic.*

**COSMO BIO CO.,LTD**

TOYO EKIMAE BLDG, 2-20, TOYO 2CHOME, KOTO-KU, TOKYO 135-0016, JAPAN  
TEL: +81-3-5632-9605 FAX: +81-3-5632-9614 e-mail: export@cosmobio.co.jp

[www.cosmobio.co.jp](http://www.cosmobio.co.jp)