

Monoclonal Antibody to Human Insulin (YO-204)

CATALOG NUMBER:	NBT-MYO-204
QUANTITY:	1mg
CLONE:	YO-204
SOURCE OF ANTIGEN:	Human Insulin
SOURCE:	BALB/C Mouse implantation ascites
FORMAT:	Purified Protein
CONCENTRATION:	10mg/ml
ISOTYPE:	IgG2a
PURIFY:	proteinA column chromatography
PURITY:	1 band (IgG1) Cellulose acetate Electrophoresis
TITER:	According to the ELISA assay, results are positive (OD 0.1) for dilutions up to 88,000 fold .
PRESENTATION:	Liquid in 10mM phosphate buffer with 0.15M NaCl , pH7.4, 0.1% Sodium Azide
STORAGE/HANDLING:	Maintain at 4 °C-8 °C (-20 °C for long term storage) in undiluted aliquots. Avoid repeated freeze/thaw cycles.

For research use only; not for use as a diagnostic.

COSMO BIO CO.,LTD

TOYO EKIMAE BLDG, 2-20, TOYO 2CHOME, KOTO-KU, TOKYO 135-0016, JAPAN
TEL: +81-3-5632-9605 FAX: +81-3-5632-9614 e-mail: export@cosmobio.co.jp

www.cosmobio.co.jp