



COSMO BIO CO., LTD.
Inspiration for Life Science

Monoclonal Antibody to Albumin (PA-062)

CATALOG NUMBER:	NBT-PA-062
QUANTITY:	1mg
CLONE:	PA-062
SOURCE OF ANTIGEN:	Albumin
SOURCE:	Rabbit searum
FORMAT:	Purified Protein
CONCENTRATION:	1mg/ml
ISOTYPE:	IgG
PURIFY:	Protein A
PURITY:	1 band (IgG1k) Cellulose acetate Electrophoresis
PRESENTATION:	Liquid in 10mM phosphate buffer with 0.15M NaCl , pH7.4, 0.1% Sodium Azide
STORAGE/HANDLING:	Maintain at 4 °C-8 °C (-20 °C for long term storage) in undiluted aliquots. Avoid repeated freeze/thaw cycles

For research use only; not for use as a diagnostic.



COSMO BIO CO., LTD.
Inspiration for Life Science

TOYO EKIMAE BLDG, 2-20, TOYO 2CHOME, KOTO-KU, TOKYO 135-0016, JAPAN

TEL: +81-3-5632-9617 FAX: +81-3-5632-9618 e-mail: export@cosmobio.co.jp

www.cosmobio.co.jp

品番. PA-062

試験成績書

抗ヒトA1bウサギ抗体

製品名: (アフィニティー精製) Lot No.

試験項目	試験結果
①特異抗体量	mg/ml (Becker 法)
②IgG濃度	1.0 mg/ml (吸光度法)
③容 量	1.0 ml
④IgG総量	1.0 mg
⑤免疫動物	ウサギ・ヤギ・ヒツジ・ニワトリ
⑥純度検定 a: 免疫電気泳動法	A1b以外の沈降線は観察されない
⑦フィブリン	外観検査では観察されない
⑧貯 法	2～8℃
⑨性 状	0.15M NaCl, 0.09% NaN ₃ を含む 10mM リン酸緩衝液(pH 7.2)に溶解

東京都国分寺市東戸倉1-15-3
株式会社日本バイオテスト研究所
T e l (042) 321-0921
F a x (042) 323-7025

コスモ・バイオ株式会社