



COSMO BIO CO., LTD.
Inspiration for Life Science

Polyclonal Antibody to Microglobulin β 2 (PG-031)

CATALOG NUMBER:	NBT-PG-031
QUANTITY:	5ml
CLONE:	PG-031
SOURCE OF ANTIGEN:	Microglobulin β 2
SOURCE:	Goat Serum
FORMAT:	Purified Protein
CONCENTRATION:	35.3 mg/ml
ISOTYPE:	IgG
PURIFY:	Protein A
PURITY:	1 band Cellulose acetate Electrophoresis
PRESENTATION:	Liquid in 10mM phosphate buffer with 0.15M NaCl , pH7.4, 0.1% Sodium Azide
STORAGE/HANDLING:	Maintain at 4 °C-8 °C (-20 °C for long term storage) in undiluted aliquots. Avoid repeated freeze/thaw cycles

For research use only; not for use as a diagnostic.



COSMO BIO CO., LTD.
Inspiration for Life Science

TOYO EKIMAE BLDG, 2-20, TOYO 2CHOME, KOTO-KU, TOKYO 135-0016, JAPAN

TEL: +81-3-5632-9617 FAX: +81-3-5632-9618 e-mail: export@cosmobio.co.jp

www.cosmobio.co.jp

品番. PG-031

試験成績書

製品名: 抗ヒト β_2 -mヤギ IgG分画 Lot. No.

試験項目	試験結果
①特異抗体量	0.2 mg/ml (Becker 法)
②蛋白濃度	mg/ml (吸光度法)
③容 量	5.0 ml
④免疫動物	ウサギ・ヤギ・セツジ・ニワトリ
⑤純度検定 a: 免疫電気泳動法	β_2 -m以外の沈降線は観察されない
⑥貯 法	2~8℃
⑦性 状	0.15M NaCl, 0.09%NaN ₃ を含む 10mMリン酸緩衝液(pH 7.2)に溶解

東京都国分寺市東戸倉1-15-3
株式会社日本バイオテスト研究所
T e l (042) 321-0921
F a x (042) 323-7025

コスモ・バイオ株式会社