



COSMO BIO CO., LTD.  
Inspiration for Life Science

## **Polyclonal Antibody to C Reactive Protein (PS-022)**

<b>CATALOG NUMBER:</b>	NBT-PS-022
<b>QUANTITY:</b>	10ml
<b>CLONE:</b>	PS-022
<b>SOURCE OF ANTIGEN:</b>	C Reactive Protein
<b>SOURCE:</b>	Rabbit
<b>FORMAT:</b>	Serum
<b>CONCENTRATION:</b>	70.0 mg/ml
<b>PRESENTATION:</b>	Liquid in 10mM phosphate buffer with 0.15M NaCl , pH7.4, 0.1% Sodium Azide
<b>STORAGE/HANDLING:</b>	Maintain at 4 °C-8 °C (-20 °C for long term storage) in undiluted aliquots. Avoid repeated freeze/thaw cycles

*For research use only; not for use as a diagnostic.*



COSMO BIO CO., LTD.  
Inspiration for Life Science

TOYO EKIMAE BLDG, 2-20, TOYO 2CHOME, KOTO-KU, TOKYO 135-0016, JAPAN  
TEL: +81-3-5632-9617 FAX: +81-3-5632-9618 e-mail: export@cosmobio.co.jp

[www.cosmobio.co.jp](http://www.cosmobio.co.jp)

品番. P S - 0 2 2

## 試験成績書

製品名：抗ヒトCRPウサギ血清      Lot No.                     

試験項目	試験結果
①特異抗体量	4. 0 mg/ml (Becker 法)
②蛋白濃度	mg/ml (屈折法)
③容 量	1 0. 0 ml
④免疫動物	ウサギ・ヤギ・ヒツジ・ニワトリ
⑤純度検定 a : 免疫電気泳動法	CRP 以外の沈降線は観察されない
⑥フィブリン	外観検査では観察されない
⑦貯 法	2 ~ 8 °C
⑧性 状	防腐剤として0.09%NaN <sub>3</sub> を含む

東京都国分寺市東戸倉1-15-3

株式会社日本バイオテスト研究所

T e l ( 0 4 2 ) 3 2 1 - 0 9 2 1

F a x ( 0 4 2 ) 3 2 3 - 7 0 2 5

**コスモ・バイオ株式会社**