



COSMO BIO CO., LTD.
Inspiration for Life Science

Polyclonal Antibody to Human Hemoglobin (PS-043)

CATALOG NUMBER:	NBT-PS-043
QUANTITY:	10ml
CLONE:	PG-043
SOURCE OF ANTIGEN:	Human Hemoglobin
SOURCE:	Sheep
FORMAT:	Serum
CONCENTRATION:	75.0 mg/ml
PURITY:	1 band Cellulose acetate Electrophoresis
PRESENTATION:	Liquid in 10mM phosphate buffer with 0.15M NaCl , pH7.4, 0.1% Sodium Azide
STORAGE/HANDLING:	Maintain at 4 °C-8 °C (-20 °C for long term storage) in undiluted aliquots. Avoid repeated freeze/thaw cycles

For research use only; not for use as a diagnostic.



COSMO BIO CO., LTD.
Inspiration for Life Science

TOYO EKIMAE BLDG, 2-20, TOYO 2CHOME, KOTO-KU, TOKYO 135-0016, JAPAN

TEL: +81-3-5632-9617 FAX: +81-3-5632-9618 e-mail: export@cosmobio.co.jp

www.cosmobio.co.jp

品番. PS-043

試験成績書

製品名: 抗ヒトHbヒツジ血清

Lot No.

試験項目	試験結果
①特異抗体量	0.7 mg/ml (Becker 法)
②蛋白濃度	mg/ml (屈折法)
③容 量	10.0 ml
④免疫動物	ウサギ・ヤギ・ヒツジ・ニワトリ
⑤純度検定 a: 免疫電気泳動法	Hb 以外の沈降線は観察されない
⑥フィブリン	外観検査では観察されない
⑦貯 法	2 ~ 8 °C
⑧性 状	防腐剤として0.09%NaN ₃ を含む

東京都国分寺市東戸倉1-15-3

株式会社日本バイオテスト研究所

Tel (042) 321-0921

Fax (042) 323-7025

コスモ・バイオ株式会社