



**POSITIVE IGNITION
WITH
ONE PUSH**



Dimensions : W 140 × D 60 × H 67 mm
Weight : 240g
Power : AA Battery

- The APT II is battery-powered, making it to be very portable. It can be used anywhere.
- The flame adjustment function of the APT II allows it to be used for a variety of purposes.
- Unique APT II spark system ensures positive ignition with one push of the button. No frustrating failure of ignition.
- The APT II was designed with compactness and safety in mind.
- A standard butane container commonly available for use in refilling cigarette lighters provides enough fuel for thousands of ignitions.
- The side window of The APT II lets you confirm quantity of gas clearly.
- The combustion pipe of bent shape protects The APT II from a trouble of wax jamming.

**Combustion
Pipe**



Check Window



COSMO BIO Co., LTD.
Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

<http://www.cosmobio.co.jp>

e-mail : export@cosmobio.co.jp

Phone : +81-3-5632-9617

FAX : +81-3-5632-9618