



COSMO BIO CO., LTD.
Inspiration for Life Science

Anti-Human ZBTB43 Polyclonal Antibody, Rabbit

Product Code: KA0414GNP

Gene: ZBTB43 (Zinc finger and BTB domain-containing protein 43, Zinc finger protein 297B, ZnF-x, Zinc finger and BTB domain-containing protein 22B)

Immunogen: GX5565 (GST-fusion protein, 155 amino acids)
QLTEHEYLPNSNSTEHDRLSTEMASQDGEEGASDSAEFHYTRPMYSK
PSIMAHKRWIHVKPERLEQACEGMDVHATYDEHQVTESINTVQTEHTV
QPSGVEEDFHIGEKKVEAEFDEQADESNYDEQVDFYGSSMEEFSGER
SDGNLIGHRQEAA

Format: Rabbit Whole Serum

Preservative : 0.05 % NaN₃

Quantity: 100 µL

Specificity: This antibody detects human ZBTB43 protein. Other species have not been tested.

Application: Immunoblotting 1:1000

Storage: -20°C

For research use only, not for diagnostic use

Manufactured by



ProteinExpress Co., Ltd. and



Kazusa DNA Res.Inst.

Kazusa DNA Research Institute.



COSMO BIO CO., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619