



COSMO BIO CO., LTD.
Inspiration for Life Science

Anti-Human BACH1 Polyclonal Antibody, Rabbit

Product Code: KB3015GNPAF

Gene: BACH1 (Transcription regulator protein BACH1, BTB and CNC homolog 1, HA2303)

Immunogen: GX5174 (GST-fusion protein, 156 amino acids)
ADQQECPRKKCFSSHQKTDLKLSSLDQRDLETDEVEEFLENKNVQTP
QCKLRRYQGNKASPPLQDSASQTYESMCLEKDAALALPSLCPKYRK
FQKAFGTDRVRTGESSVKDIHASVQPNRSENECLGGVPECRDLQVM
LKCDESKLAMEPEE

Format: Affinity-purified rabbit IgG solubilized in PBS containing with 40% glycerol and 0.02% of NaN₃.

Quantity: 50 µg (250 µL)

Specificity: This antibody detects human BACH1 protein. Other species have not been tested.

Application: Immunoblotting 1:1000

Storage: -20°C

For research use only, not for diagnostic use

Manufactured by



ProteinExpress Co., Ltd. and



Kazusa DNA Research Institute.

Kazusa DNA Res.Inst.



COSMO BIO CO., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619