



COSMO BIO CO., LTD.
Inspiration for Life Science

Anti-Human LMO4 Polyclonal Antibody, Rabbit

Product Code: KB3107GNP

Gene: LMO4 (LIM domain transcription factor LMO4, LIM domain only protein 4, LMO-4, Breast tumor autoantigen)

Immunogen: GX5086 (GST-fusion protein, 163 amino acids)
MVNPGSSSQPPPVTAGSLSWKRCAGCGGKIADRFLLYAMDSYWHSRC
LKSCCQAQLGDIGTSCYTKSGMILCRNDYIRLFGNSGACSACGQSIPA
SELMRAQGNVYHLKCTCSTCRNRLVPGDRFHYINGSLFCEHDRPTA
LINGHLNSLQSNPLLPDQK

Format: Rabbit Whole Serum

Preservative : 0.05 % NaN₃

Quantity: 100 µL

Specificity: This antibody detects human LMO4 protein. Other species have not been tested.

Application: Immunoblotting 1:1000

Storage: -20°C

For research use only, not for diagnostic use

Manufactured by



ProteinExpress Co., Ltd. and



Kazusa DNA Res.Inst.

Kazusa DNA Research Institute.



COSMO BIO CO., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619