



COSMO BIO CO., LTD.
Inspiration for Life Science

Anti-Human CDK7 Polyclonal Antibody, Rabbit

Product Code: KB3201GNPAF

Gene: CDK7 (Cell division protein kinase 7, CDK-activating kinase, CAK, CAK1, TFIIH basal transcription factor complex kinase subunit, 39 kDa protein kinase, P39 Mo15, STK1)

Immunogen: GX5177 (GST-fusion protein, 173 amino acids)
MALDVKSRAKRYEKLDLFLGEGQFATVYKARDKNTNQIVAIIKKIKLGHRS
EAKDGINRTALREIKLLQELSHPNIIIGLLDAFGHKSNI SLVDFMETDLEVI
IKDNSLVLT PSHIKAYMLMTLQGLEYLHQHWILHRDLKPNNLLLDENGVL
KLADFGLA KSFGSPN RAYTHQV

Format: Affinity-purified rabbit IgG solubilized in PBS containing with 40% glycerol and 0.02% of NaN₃.

Quantity: 50 µg (250 µL)

Specificity: This antibody detects human CDK7 protein. Other species have not been tested.

Application: Immunoblotting 1:1000

Storage: -20°C

For research use only, not for diagnostic use

Manufactured by



ProteinExpress Co., Ltd. and



Kazusa DNA Res.Inst.

Kazusa DNA Research Institute.



COSMO BIO CO., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619