



COSMO BIO CO., LTD.
Inspiration for Life Science

Anti-Human FEM1A Polyclonal Antibody, Rabbit

Product Code: KB3219GNPAF

Gene: FEM1A (Protein fem-1 homolog A, FEM1-alpha, FEM1a, Prostaglandin E receptor 4-associated protein)

Immunogen: GX5237 (GST-fusion protein, 157 amino acids)
APCCSSSPPEEPLNGESYESCPTSREAAVEALELLGATYVDKKRDLLG
ALKHWRRAMELRHQGGEYLPKPEPPQLVLAYDYSREVNTTEELEALIT
DPDEMRMQALLIRERILGPSHPDTSYYIRYRGAVYADSGNFERCIRLWK
YALDMQQSNLEP

Format: Affinity-purified rabbit IgG solubilized in PBS containing with 40% glycerol and 0.02% of NaN₃.

Quantity: 50 µg (250 µL)

Specificity: This antibody detects human FEM1A protein. Other species have not been tested.

Application: Immunoblotting 1:1000

Storage: -20°C

For research use only, not for diagnostic use

Manufactured by



ProteinExpress Co., Ltd. and



Kazusa DNA Res.Inst.

Kazusa DNA Research Institute.



COSMO BIO CO., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619