



COSMO BIO CO., LTD.
Inspiration for Life Science

Anti-Human CBX5 Polyclonal Antibody, Rabbit

Product Code: KB3469GNPAF

Gene: CBX5 (Chromobox protein homolog 5, Heterochromatin protein 1 homolog alpha, HP1 alpha, Antigen p25)

Immunogen: GX5171 (GST-fusion protein, 161 amino acids)
YVVEKVLDRRVVKGQVEYLLKWKGFSSEHNTWEPEKNLDCPELISEFM
KKYKKMKEGENNKPRESSESNNRKSNSNSADDIKSKKKREQSNDIAR
GFERGLEPEKIIGATDSCGDLMLMKWKDTDEADLVLAKEANVKCPQIV
IAFYEEERLTWHAYPED

Format: Affinity-purified rabbit IgG solubilized in PBS containing with 40% glycerol and 0.02% of NaN₃.

Quantity: 50 µg (250 µL)

Specificity: This antibody detects human CBX5 protein. Other species have not been tested.

Application: Immunoblotting 1:1000

Storage: -20°C

For research use only, not for diagnostic use

Manufactured by



ProteinExpress Co., Ltd. and



Kazusa DNA Res.Inst.

Kazusa DNA Research Institute.



COSMO BIO CO., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619