



COSMO BIO CO., LTD.
Inspiration for Life Science

Anti-Human MEF2A Polyclonal Antibody, Rabbit

Product Code: KB3473GNP

Gene: MEF2A (Myocyte-specific enhancer factor 2A, Serum response factor-like protein 1)

Immunogen: GX5088 (GST-fusion protein, 165 amino acids)
GANSLGKVMPTKSPPPPGGGNLGMNSRKPDLRVVIPPSSKGMMPPPLN
TQRISSQATQPLATPVVSVTTPSLPPQGLVYSAMPTAYNTDYSLTSD
LSALQGFNSPGMLSLGQVSAWQQHHLGQAALSSLVAGGQLSQGSNLS
INTNQNISIKSEPIPPRDRMT

Format: Rabbit Whole Serum

Preservative : 0.05 % NaN₃

Quantity: 100 µL

Specificity: This antibody detects human MEF2A protein. Other species have not been tested.

Application: Immunoblotting 1:1000

Storage: -20°C

For research use only, not for diagnostic use

Manufactured by



ProteinExpress Co., Ltd. and



Kazusa DNA Res.Inst.

Kazusa DNA Research Institute.



COSMO BIO CO., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619