



COSMO BIO CO., LTD.
Inspiration for Life Science

Anti-Human ETV1 Polyclonal Antibody, Rabbit

Product Code: KB3955GNP

Gene: ETV1 (ETS translocation variant 1, Protein ER81)

Immunogen: GX5330 (GST-fusion protein, 172 amino acids)
DQKPQVGMRPSNPPTPSSTPV SPLHHASPNSTHTPKPDRAFP AHLPP
SQSIPDSSYPMDHRFRRLSEPCNSFPPLPTMPREGRPMYQRQMSEP
NIPFPPQGFKQEYHDPVYEHTMVGSAASQSFPPPLMIKQEPRDFAYD
SGCMFEKGPRQFYDDTCVVPEKFDGDIKQE

Format: Rabbit Whole Serum

Preservative : 0.05 % NaN₃

Quantity: 100 µL

Specificity: This antibody detects human ETV1 protein. Other species have not been tested.

Application: Immunoblotting 1:1000

Storage: -20°C

For research use only, not for diagnostic use

Manufactured by



ProteinExpress Co., Ltd. and



Kazusa DNA Res.Inst.

Kazusa DNA Research Institute.



COSMO BIO CO., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619