



COSMO BIO CO., LTD.
Inspiration for Life Science

Anti-Human PBX1 Polyclonal Antibody, Rabbit

Product Code: KB4019GNPAF

Gene: PBX1 (Pre-B-cell leukemia transcription factor 1, Homeobox protein PBX1, Homeobox protein PRL)

Immunogen: GX5227 (GST-fusion protein, 174 amino acids)
MDEQPRLMHSAGVGMAGHPGLSQHLQDGAGGTEGEGGRKQDIGDI
LQQIMTITDQSLDEAQARKHALNCHRMKPALFNVLCIEKEKTVLSIRGA
QEEPTDPQLMRLDNMLLAEGVAGPEKGGGSAAAAAASGGAGS
DNSVEHSDYRAKLSQIRQIYHTELEKYEACNE

Format: Affinity-purified rabbit IgG solubilized in PBS containing with 40% glycerol and 0.02% of NaN₃.

Quantity: 50 µg (250 µL)

Specificity: This antibody detects human PBX1 protein. Other species have not been tested.

Application: Immunoblotting 1:1000

Storage: -20°C

For research use only, not for diagnostic use

Manufactured by



ProteinExpress Co., Ltd. and



Kazusa DNA Res.Inst.

Kazusa DNA Research Institute.



COSMO BIO CO., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619