



COSMO BIO CO., LTD.
Inspiration for Life Science

Anti-Human SMARCC2 Polyclonal Antibody, Rabbit

Product Code: KB4171GNPAF

Gene: SMARCC2 (SWI/SNF complex subunit SMARCC2, SWI/SNF-related matrix-associated actin-dependent regulator of chromatin subfamily C member 2, SWI/SNF complex 170 kDa subunit, BRG1-associated factor 170)

Immunogen: GX5112 (GST-fusion protein, 170 amino acids)
MAVRKKDGGPNVKYYEAADTVTQFDNVRLWLGNKYKKYIQAEPTNK
SLSSLVVQLLQFQEEVFGKHVSNAPLTKLPIKCFLDFKAGGSLCHILAAA
YKFKSDQGWRRYDFQNPSRMDRNVEMFMTIEKSLVQNNCLSRPNIFL
CPEIEPKLLGKLKDIIKRHQGTVED

Format: Affinity-purified rabbit IgG solubilized in PBS containing with 40% glycerol and 0.02% of NaN₃.

Quantity: 50 µg

Specificity: This antibody detects human SMARCC2 protein. Other species have not been tested.

Application: Immunoblotting 1:1000

Storage: -20°C

For research use only, not for diagnostic use

Manufactured by



ProteinExpress Co., Ltd. and



Kazusa DNA Res.Inst.

Kazusa DNA Research Institute.



COSMO BIO CO., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619