



COSMO BIO CO., LTD.
Inspiration for Life Science

Anti-Human MLLT6 Polyclonal Antibody, Rabbit

Product Code: KB4336GNP

Gene: MLLT6 (myeloid/lymphoid or mixed-lineage leukemia (trithorax homolog, Drosophila); translocated to, 6, Protein AF-17)

Immunogen: GX5064 (GST-fusion protein, 173 amino acids)
EEDKYSKPTAPAPSAPPSPSAPEPPKADLFEQKVVFSGFGPIMRFSTTT
SSSGRARAPSPGDYKSPHVTGSGASAGTHKMPALSATPVPADETPE
TGLKEKKHKASKRSRHGPRPKGSRNKEGTGGPAAPSLPSAQLAGFT
ATAASPFSGGSLVSSGLGGLSSRTFGPSGS

Format: Rabbit Whole Serum

Preservative : 0.05 % NaN₃

Quantity: 100 µL

Specificity: This antibody detects human MLLT6 protein. Other species have not been tested.

Application: Immunoblotting 1:1000

Storage: -20°C

For research use only, not for diagnostic use

Manufactured by



ProteinExpress Co., Ltd. and



Kazusa DNA Res.Inst.

Kazusa DNA Research Institute.



COSMO BIO CO., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619