



COSMO BIO CO., LTD.
Inspiration for Life Science

Anti-Human PML Polyclonal Antibody, Rabbit

Product Code: KB4394GNP

Gene: PML (Probable transcription factor PML, Tripartite motif-containing protein 19, RING finger protein 71)

Immunogen: GX5060 (GST-fusion protein, 173 amino acids)
ARLARSSPEQPRPSTSKAVSPPHLDGPPSPRSPVIGSEVFLPNSNHVA
SGAGEAEERVVVISSSESDSAENSSSRELDSSSESSDLQLEGPSTLR
VLDENLADPQAEDRPLVFFDLKIDNETQKISQLAAVNRESKFRVVIQPEA
FFSIYSKAVSLEVGLQHFLSFLSSMRR

Format: Rabbit Whole Serum

Preservative : 0.05 % NaN₃

Quantity: 100 µL

Specificity: This antibody detects human PML protein. Other species have not been tested.

Application: Immunoblotting 1:1000

Storage: -20°C

For research use only, not for diagnostic use

Manufactured by



ProteinExpress Co., Ltd. and



Kazusa DNA Res.Inst.

Kazusa DNA Research Institute.



COSMO BIO CO., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619