



COSMO BIO CO., LTD.
Inspiration for Life Science

Anti-Human ZNF238 Polyclonal Antibody, Rabbit

Product Code: KB4680GNP

Gene: ZNF238 (Zinc finger protein 238, Transcriptional repressor RP58, 58kDa repressor protein, Zinc finger protein C2H2-171, Translin-associated zinc finger protein 1, TAZ-1, Zinc finger and BTB domain-containing protein 18)

Immunogen: GX5181 (GST-fusion protein, 150 amino acids)
DLAAEPGNMWMRLPSDSAGIPQAGGEAEPHATAAGKTVASPCSSSTESLSQRSVTSVRD
SADVDCVLDLSVKSSLSGVENLNSSYFSSQDVLRSNLVQVKVEKEASCDESVDGTNDY
DMEHSTVKESVSTNNRVQYEPAPHLAPLREDSVLR
ROUGE Database (Mouse KIAA/FLJ cDNA Database)
[GX5181](#)
The mammalian expression vector of the full length human protein corresponding to the antibody is available!
[Human cDNA/ORF Clones in Kazusa Collection \(KDRI\)](#)

Format: Rabbit Whole Serum

Preservative : 0.05 % NaN₃


Quantity: 100 µL

Specificity: This antibody detects human ZNF238 protein. Other species have not been tested

Application: Western blotting (1 : 1,000), Other applications have not been tested.

Storage: Store at -20°C. Avoid freeze-thaw cycles

Anti-Human ZNF238



For research use only, not for diagnostic use

Manufactured by



ProteinExpress Co., Ltd. and



Kazusa DNA Research Institute.



COSMO BIO CO., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.com>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619