



COSMO BIO CO., LTD.
Inspiration for Life Science

Anti-Human STAT6 Polyclonal Antibody, Rabbit

Product Code: KB4980GNP

Gene: STAT6 (Signal transducer and activator of transcription 6, IL-4 Stat)

Immunogen: GX5076 (GST-fusion protein, 147 amino acids)
RDLAQLKNLYPKPKDEAFRSHYKPEQMGKDGRGYVPATIKMTVERD
QPLPTPELQMPTMVPSYDLGMAPDSSMSMQLGPDMVPQVYPPHSHSI
PPYQGLSPEESVNVLSAFQEPHLQMPPSLGQMSLPFDQPHPQGLLPC
QPQEHA

Format: Rabbit Whole Serum

Preservative : 0.05 % NaN₃

Quantity: 100 µL

Specificity: This antibody detects human STAT6 protein. Other species have not been tested.

Application: Immunoblotting 1:1000

Storage: -20°C

For research use only, not for diagnostic use

Manufactured by



ProteinExpress Co., Ltd. and



Kazusa DNA Research Institute.

Kazusa DNA Res.Inst.



COSMO BIO CO., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619