



COSMO BIO CO., LTD.
Inspiration for Life Science

Anti-Human STAT3 Polyclonal Antibody, Rabbit

Product Code: KB5090GNPAF

Gene: STAT3 (Signal transducer and activator of transcription 3, Acute-phase response factor)

Immunogen: GX5069 (GST-fusion protein, 167 amino acids)
AILSTKPPGTFLLRFSESSKEGGVTFTWVEKDISGKTQIQSVEPYTKQQ
LNNMSFAEIIMGYKIMDATNILVSPLYLYPDIPKEEAFGKYCRPESQEHP
EADPGAAPYLKTKFICVPTTCSNTIDLPMSPRTLDSLMLQFGNNGEGAE
PSAGGQFESLTFDMELTS

Format: Affinity-purified rabbit IgG solubilized in PBS containing with 40% glycerol and 0.02% of NaN₃.

Quantity: 50 µg (250 µL)

Specificity: This antibody detects human STAT3 protein. Other species have not been tested.

Application: Immunoblotting 1:1000

Storage: -20°C

For research use only, not for diagnostic use

Manufactured by



ProteinExpress Co., Ltd. and



Kazusa DNA Res.Inst.

Kazusa DNA Research Institute.



COSMO BIO CO., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619